


Immunodox®-Maleimide (PEGylated) (Post-insertion)

DESCRIPTION

During the past five decades, various types of chemistries have been used for conjugation of molecules such as antibodies to the surface of the liposomes. In general, the conjugation can be achieved through the N-terminus, the C-terminus or the available sulfur (*e.g.* Fab' fraction or thiolated Ab). Not all chemistries have the same yield and efficiency of conjugation and often reproducing biocompatible batches can be a challenge. Coupling of sulfhydryl groups with maleimide groups has been the most widely used conjugation of antibodies to liposomes. Different lipids which are offered for thioether conjugation contain maleimide, aromatic maleimides such as N-[4-(p-maleimidophenyl)-butyryl] (MPB) or 4-(N-maleimidomethyl)cyclohexane-1-carboxylate (MCC) group. The maleimide function group of MCC which contains an aliphatic cyclohexane ring is more stable toward hydrolysis in aqueous reaction environments rather than the aromatic phenyl group of MPB. MPB and MCC lipids are non-PEGylated lipids and they have separate kits and protocols than PEGylated maleimide lipids.

One of the major problems of using maleimide chemistry for conjugation is the rapid hydrolysis of maleimide lipid. The rate of hydrolysis is much faster in alkaline pH and therefore controlling the pH throughout the entire process is necessary and it is recommended to use the pH of 7. Due to the hydrolysis of maleimide group, our kits are designed for post-insertion of ligand conjugated maleimide lipid into the preformed liposomes. After post conjugation the liposomes must be used right away because hydrolysis may occur after sulfhydryl coupling to the maleimide as well. Another problem is the reactivity and oxygen sensitivity of sulfhydryl group on thiolated antibody or Fab' fragment. Due to that the conjugation reaction should be done under argon or nitrogen using inflatable polyethylene glove bag chambers.

Thiolation which is adapted to the modification of all the antibody functional groups, is relatively clean, fast, and efficient. However, different antibodies may be more sensitive to some procedures than others. Therefore, it is recommended to select the chemistry and site of modification depending on what procedures are compatible with the antibody.



Conjugation reaction between maleimide-activated DSPE-PEG lipid with the sulfhydryl group of the ligand. The micelles formed from lipid conjugated ligand and non-reactive PEG lipids are mixed together and the PEGylated lipids are post-inserted into the liposomes to form PEGylated ligand surface conjugated liposomes.

FORMULATION INFORMATION

Immunodox®-Maleimide (PEGylated) (Post-insertion)

Post-insertion Kit (3 Vials)	Specification
Vial 1	Preformed Liposomes composed of HSPC and Cholesterol (60:40 molar ratio)
Vial 2	DSPE-PEG(2000)-Maleimide lipid (reactive PEGylated lipid) in powder form
Vial 3	DSPE-PEG(2000) lipid (non-reactive PEGylated lipid) in powder form

Lipid Composition for Vial 1*	Concentration (mg/ml)	Concentration (mM)	Molar Ratio Percentage
Hydrogenated Soy PC	11.5	14.66	60
Cholesterol	3.83	9.9	40
Total	15.33	24.56	100


* For the 5-ml kit, the volume of vial 1 is 4 ml. 1 ml of micelle solution that are formed using vials 2 and 3 will be added to this vial to make the final volume of 5 ml in the final product. For the 2-ml kit, the volume of vial 1 is 1.6 ml. 0.4 ml of micelle solution that is formed using vials 2 and 3 will be added to this vial to make the final volume of 2 ml in the final product.

Buffers, Liposome Size and Encapsulated Drug Concentration	Specification
Inside Buffer	Ammonium Sulfate
Outside Buffer	Phosphate Buffered Saline
pH	7.4
Liposome Size	100 nm
Encapsulated Doxorubicin	2 mg/ml (3.45 mM)

Vial 2*	Specification
---------	---------------

DSPE-PEG(2000)-Maleimide Lipid

This vial contains reactive DSPE-PEG(2000)-Maleimide lipid in powder form. This lipid is conjugated to a reactive protein, peptide or ligand containing sulfhydryl and then mixed with non-reactive DSPE-PEG(2000) lipid in aqueous solution to form micelles. The PEGylated lipid micelles are incubated with preformed liposomes in vial 1 and PEG lipids will post-insert themselves into the liposomes.




* The amount of the powdered PEG(2000)-Maleimide lipid for the 2-ml kit is 1.34 mg and for the 5-ml kit is 3.34 mg.

Vial 3*
Specification

 DSPE-PEG(2000)
Lipid

This vial contains non-reactive DSPE-PEG(2000) lipid in powder form. This lipid in mixed with DSPE-PEG(2000)-Maleimide lipid which is already conjugated to a ligand (protein, peptide, *etc.*) in aqueous solution to form micelles. The PEGylated lipid micelles are incubated with preformed liposomes in vial 1 and PEG lipids will post-insert themselves into the liposomes.



* The amount of the powdered PEG(2000)-DSPE lipid for the 2-ml kit is 5 mg and for the 5-ml kit is 12.5 mg.

CONJUGATION PROTOCOL (POST-INSERTION)

Materials and Equipment

The 3-vial post-insertion kit contains preformed liposomes (vial 1), DSPE-PEG(2000)-Maleimide lipid in powder form (vial 2) and non-reactive PEGylated lipid in powder form (vial 3). In order to use the post-insertion kit, you will need:

1. Two small 10-ml round bottom flasks or two small glass vials.
2. A rotary evaporator. We understand that many labs might not have a rotovap. Alternatively, you can use a nitrogen tank connected to a thin hose for creating a stream of nitrogen flow to dry the lipid and make a thin film.

3. A small amount of a solvent such as chloroform or methylene chloride (you will only need a few milliliters).
4. Phosphate buffered saline (PBS). pH should be adjusted to 7.
5. [2-mercaptoethanol](#).
6. [Aldrich®-Atmosbag](#) connected to a nitrogen tank. Due to oxygen sensitivity of the reaction, the coupling reaction should be done in an oxygen-free environment.
7. [Float-A-Lyzer®](#) with a proper MWCO that easily allows the cleanup of your liposome conjugated ligand from free and non-conjugated protein/peptide/ligand. You need to make sure that the MWCO is below 1,000,000 dalton. At 1,000,000 dalton, the pore size on the dialysis membrane gets close to 100 nm and therefore, your liposomes can be dialyzed out. You cannot use dialysis cassettes blindly. Please understand the technique before using either spin columns or dialysis cassettes. If you do not use the correct MWCO, you can lose your entire prep. For this protocol, we recommend MWCO of 300,000 dalton.
8. A Sonicator. It is better to have a bath sonicator. If you do not, that is fine, and you still can follow the protocol. You may also use a vortex instead of the sonicator for agitation of the solution as well.

Preparation Method

1. The post-insertion kits come in two sizes; 2 ml and 5 ml. For the 2-ml kit, dissolve the content of vial 3 (non-reactive PEGylated lipid) in 100 μ l of chloroform or methylene chloride. For the 5-ml kit, the content of vial 3 should be dissolved in 250 μ l of chloroform or methylene chloride. Transfer the solution to a 10-ml round bottom flask. Dry the chloroform using a rotary evaporator or under a stream of nitrogen to make a dried lipid film.
2. For the 2-ml kit, add 100 μ l of PBS buffer to the dried lipid film. For the 5-ml kit, the amount of the added buffer is 250 μ l. It is preferred to sonicate the hydrated lipid film using a bath sonicator and sonicate the micelle solution for 5 minutes. If you do not have

a bath sonicator, then hydrate the dried lipid film with PBS for at least 1 hour and constantly rotate the solution in the round bottom flask using a rotavap (not connected to vacuum) or by hand to make sure that all the dried lipid on the wall of the round bottom flask will go to the solution and form micelles. Alternatively, you can use a vortex to agitate the solution. The goal is to have all the dried lipid on the wall of the round bottom glass to go to the micelle solution. Cover the mouth of the round bottom flask with parafilm. Refrigerate the micelle solution of non-reactive PEG lipids until it is ready to be mixed with micelles formed in the step 5.

3. The 2-ml kit contains 1.30 mg (0.22 μmol) of reactive DSPE-PEG(2000)-Maleimide lipid (vial 2). The 5-ml kit contains 3.25 mg (0.55 μmol) of reactive DSPE-PEG(2000)-Maleimide lipid (vial 2). For the 5-ml kit size, the content of vial 2 (DSPE-PEG(2000)-Maleimide lipid) should be dissolved in 250 μl of chloroform or methylene chloride. Transfer the solution to a 10-ml round bottom flask. Dry the chloroform using a rotary evaporator or under a stream of nitrogen to make a dried lipid film.
4. Dried DSPE-PEG-Maleimide film is hydrated with PBS buffer to form a micellar lipid solution. If you are using the 2-ml post-insertion kit, then hydrate the 1.30 mg of dried DSPE-PEG-Maleimide lipid film in 100 μl of buffer, and if you are using the 5-ml post-insertion kit, then hydrate the 3.25 mg of dried DSPE-PEG-Maleimide lipid film in 250 ml of buffer.
5. Incubate the micellar lipid solution with the antibody, protein or peptide at 3:1 molar ratio or lipid to protein. Allow the reaction to proceed in phosphate buffer under the nitrogen (inert gas) chamber for 8 hours at room temperature with moderate stirring. The concentration of antibody, peptide or protein that is added to micellar solution is depend on the solubility of your molecule. It is recommended to use a fairly concentrated solution. For example, use a volume around 100 μl of antibody, peptide or protein for the 2-ml kit and around 250 μl of antibody, peptide or protein for the 5-ml kit.
6. The excess maleimide groups were capped by reaction with [2-mercaptoethanol](#). The reaction is quenched with 2 mM 2-mercaptoethanol for 30 min.

7. The micelles obtained from the steps 2 and 5 are mixed. Total volume of the 2 mixed micelles for the 2-ml kit is 300 μ l and for the 5-ml kit is 750 μ l. Incubate the mixed micelles with preformed liposomes (vial 1) at 60 °C for 30 min.
8. Remove non-conjugated antibody, protein, peptide or ligand by dialysis. We prefer dialysis to size exclusion columns. Dialysis is a much slower process but there will be minimum loss of immunoliposomes after the prep is cleaned from non-conjugated protein/peptide/ligand. Spin columns are much faster, but you can easily lose over 50% of the liposomes on the spin column. We recommend using [Float-A-Lyzer®](#) dialysis cassette from Spectrum Labs. You need to choose a cassette with proper MWCO depending on the MW of your protein, ligand, antibody or antibody fragment. In this case, we recommend using a dialysis cassette with MWCO of 300,000 dalton. **NOTE:** If you decide to use a dialysis cassette, you need to make sure that the MWCO is below 1,000,000 dalton. At 1,000,000 dalton, the pore size on the dialysis membrane gets close to 100 nm and therefore your liposomes can be dialyzed out. You cannot use dialysis cassettes and spin columns blindly. They come in various sizes, and you need to choose the correct size wisely. Dialyze the immunoliposome solution in 1 liter of PBS at pH 7 for 8 hours. Change the dialysis buffer with a fresh 1 liter of PBS and let it dialyze for another 8 hours. After this step, your cleaned up immunoliposome is ready to be used.


QUANTIFICATION OF REACTIVE SULFHYDRYL IN ANTIBODIES OR LIGANDS (ELLMAN'S ASSAY)

The yield of conjugation is the most important factor in formulating immunoliposomes. Many scientists simply assume that their thiolated antibody or the Fab' fraction contains reactive sulfhydryl for conjugation to maleimide lipid without further assay. Disulfide bridge can form very easily so it is very important to quantify the available reactive sulfhydryl in your antibody or ligand solution before performing the conjugation reaction with maleimide liposomes.

Ellman's assay is a widely used assay for determining the amount of free sulfhydryl. You can follow the step by step protocol [here](#).

TECHNICAL NOTES

- Doxorubicin is a fluorescent molecule with λ_{ex} 470 nm and λ_{em} 585 nm. If you are using a fluorescent tag on your antibody or ligand, you need to make sure that they will not interfere with each other.
- After conjugation reactions, liposomes containing excess maleimide or thiol groups may exhibit undesirable qualities, such as aggregation, reactions *in vitro* and *in vivo*, and immunogenicity. These reactive moieties can be quenched with reagents containing iodo-, maleimide, or sulfhydryl groups where appropriate. This is likely to be a particularly serious problem for thiolated liposomes. Therefore, it is recommended that the antibody to be thiolated to generate the appropriate reactive entities for the final conjugation reaction.
- In order to prevent oxidation of sulfhydryl on antibody and formation of disulfide bridge, the coupling reaction must be performed under an inert atmosphere such as argon or nitrogen. To set up an inert gas chamber we recommend using [Aldrich®-Atmosbag](#) which is a flexible, inflatable polyethylene chamber with built-in gloves which is a portable and inexpensive alternative to laboratory glove box.
- Maleimide group on lipid is highly sensitive of alkaline pH and it will hydrolyze rapidly at higher pH. Experimental investigations have been shown that in alkaline condition (pH > 7.5), maleimide and its derivatives are hydrolyzed to a non-reactive maleamic acid (see the figure below). This instability should be considered in any quantitative procedures, such as coupling with sulfhydryl groups. Therefore, it is very important to make sure that the pH of the reaction with stay between 6.5 and 7 during the entire process.
- Liposomes should be kept at 4 °C and **NEVER** be frozen.



Dissociation curve for the equilibrium between maleimide and its anion.

(adapted from reference 2)

APPEARANCE

Immunodox®-Maleimide is a red translucent liquid made of nano size unilamellar liposomes. Usually due to the small size of liposomes no settling will occur in the bottom of the vial. The liposomes are packaged in an amber vial.

STORAGE AND SHELF LIFE

Storage

Immunodox® products should always be stored at in the dark at 4 °C, except when brought to room temperature for brief periods prior to animal dosing. DO NOT FREEZE. If the suspension is frozen, the encapsulated drug can be released from the liposomes thus limiting its effectiveness. In addition, the size of the liposomes will also change upon freezing and thawing.

Shelf Life

Immunodox®-Maleimide is made on daily basis. The batch that is shipped is manufactured on the same day. It is advised to use the products within 4 months of the manufacturing date.

REFERENCES AND BACKGROUND READING

1. *Matsui, S., and H. Aida. "Hydrolysis of some N-alkylmaleimides." Journal of the Chemical Society, Perkin Transactions 2 12 (1978): 1277-1280.*
2. *Barradas, Remigio Germano, Stephen Fletcher, and John Douglas Porter. "The hydrolysis of maleimide in alkaline solution." Canadian Journal of Chemistry 54.9 (1976): 1400-1404.*
3. *Gregory, John D. "The stability of N-ethylmaleimide and its reaction with sulfhydryl groups." Journal of the American Chemical Society 77.14 (1955): 3922-3923.*
4. *Nassander UK, Steerenberg PA, De Jong WH, Van Overveld WO, Te Boekhorst CM, Poels LG, Jap PH, Storm G. Design of immunoliposomes directed against human ovarian carcinoma. Biochimica et Biophysica Acta (BBA)-Biomembranes. 1995 Apr 12;1235(1):126-39.*
5. *Derksen JT, Morselt HW, Scherphof GL. Uptake and processing of immunoglobulin-coated liposomes by subpopulations of rat liver macrophages. Biochimica et Biophysica Acta (BBA)-Molecular Cell Research. 1988 Sep 16;971(2):127-36.*

6. *Derksen JT, Scherphof GL. An improved method for the covalent coupling of proteins to liposomes. Biochimica et Biophysica Acta (BBA)-Biomembranes. 1985 Mar 28;814(1):151-5.*
7. *Nässander UK, Steerenberg PA, Poppe H, Storm G, Poels LG, De Jong WH, Crommelin DJ. In vivo targeting of OV-TL 3 immunoliposomes to ascitic ovarian carcinoma cells (OVCAR-3) in athymic nude mice. Cancer research. 1992 Feb 1;52(3):646-53.*
8. *Park JW, Hong K, Carter P, Asgari H, Guo LY, Keller GA, Wirth C, Shalaby R, Kotts C, Wood WI. Development of anti-p185HER2 immunoliposomes for cancer therapy. Proceedings of the National Academy of Sciences. 1995 Feb 28;92(5):1327-31.*
9. *Koning GA, Morselt HW, Velinova MJ, Donga J, Gorter A, Allen TM, Zalipsky S, Kamps JA, Scherphof GL. Selective transfer of a lipophilic prodrug of 5-fluorodeoxyuridine from immunoliposomes to colon cancer cells. Biochimica et Biophysica Acta (BBA)-Biomembranes. 1999 Aug 20;1420(1):153-67.*
10. *Vingerhoeds MH, Steerenberg PA, Hendriks JJ, Dekker LC, Van Hoesel QG, Crommelin DJ, Storm G. Immunoliposome-mediated targeting of doxorubicin to human ovarian carcinoma in vitro and in vivo. British journal of cancer. 1996 Oct 1;74(7):1023-9.*
11. *Kirpotin D, Park JW, Hong K, Zalipsky S, Li WL, Carter P, Benz CC, Papahadjopoulos D. Sterically stabilized anti-HER2 immunoliposomes: design and targeting to human breast cancer cells in vitro. Biochemistry. 1997 Jan 7;36(1):66-75.*
12. *Garnier B, Bouter A, Gounou C, Petry KG, Brisson AR. Annexin A5-functionalized liposomes for targeting phosphatidylserine-exposing membranes. Bioconjugate chemistry. 2009 Oct 19;20(11):2114-22.*
13. *Mattson G, Conklin E, Desai S, Nielander G, Savage MD, Morgensen S. A practical approach to crosslinking. Molecular biology reports. 1993 Apr 1;17(3):167-83.*
14. *Smyth DG, Blumenfeld OO, Konigsberg W. Reactions of N-ethylmaleimide with peptides and amino acids. Biochemical Journal. 1964 Jun;91(3):589-95.*

15. Harokopakis E, Childers NK, Michalek SM, Zhang SS, Tomasi M. Conjugation of cholera toxin or its B subunit to liposomes for targeted delivery of antigens. Journal of immunological methods. 1995 Sep 11;185(1):31-42.
16. Gradauer K, Vonach C, Leitinger G, Kolb D, Fröhlich E, Roblegg E, Bernkop-Schnürch A, Prassl R. Chemical coupling of thiolated chitosan to preformed liposomes improves mucoadhesive properties. International journal of nanomedicine. 2012;7:2523-34.
17. Gradauer K, Barthelmes J, Vonach C, Almer G, Mangge H, Teubl B, Roblegg E, Dünnhaupt S, Fröhlich E, Bernkop-Schnürch A, Prassl R. Liposomes coated with thiolated chitosan enhance oral peptide delivery to rats. Journal of controlled release. 2013 Dec 28;172(3):872-8.
18. Moreira JN, Ishida T, Gaspar R, Allen TM. Use of the post-insertion technique to insert peptide ligands into pre-formed stealth liposomes with retention of binding activity and cytotoxicity. Pharmaceutical research. 2002 Mar 1;19(3):265-9.