

## Five-Color Mannosylated Fluoroliposome® Kit

### DESCRIPTION

---


Mannose receptor targeting by mannosylated liposomes has been demonstrated for a variety of mannosylated lipid conjugates in a variety of liposome morphologies and compositions in several different *in vitro* and *in vivo* models. A very large number of publications is about using a hydrophobic derivative of mannose (4-aminophenyl  $\alpha$ -D-mannopyranoside) rather than using a mannosylated lipid in clodronate liposomes. This is mainly due to the high cost and complexity of synthesizing and conjugating mannose to lipid. 4-aminophenyl  $\alpha$ -D-mannopyranoside is commercially available and far less expensive than synthesizing mannose conjugated lipid.

Why mannose? Mannose is one of the carbohydrate components of many bacterial and viral cell surfaces; therefore, the ever-efficient, highly redundant immune system has evolved multiple mechanisms for identifying pathogens based on mannose recognition. The animal and plant kingdoms likewise utilize carbohydrate recognition signaling mechanisms including mannose residues. Many publications evaluate other carbohydrates as targeting mechanisms for various cell types, however mannose targeting to phagocytes appears to be one of the more specific mechanisms identified to date. Mammalian cell surface identification molecules based on mannose binding, such as the ICAM family of leukocyte adhesion molecules, target the SIGN family of mannose receptors to accomplish self-recognition *in vivo*.


A well-known and cited study by Eto *et al.* [1] demonstrates that liposomes containing aminophenyl mannoside were most efficiently incorporated into the mouse brain across the blood brain barrier. The radiolabeled liposomes bearing aminophenyl  $\alpha$ -D-mannopyranoside were maximally incorporated into the mouse brain after 48 hours, whereas in the spleen and liver, these radioactivities were maximum after 12 hours. The studies also showed that liposomes were most incorporated was glial cells rather than neuronal cell. The subcellular fractionation study indicates that mannose labeled liposomes are incorporated into lysosomes rich fraction both in liver and brain.

There are five mannosylated fluorescent control liposome products (m-Fluoroliposome®) for m-Clodrosome® (mannosylated clodronate liposomes). All five mannosylated fluorescent liposomes incorporate a lipophilic dye inside their membranes. They are insoluble in water; however, their fluorescence is easily detected when incorporated into membranes. DiI, DiO, DiD, DiR and DiA cover a wide range of excitation and emission wavelengths from 300s to 900s. DiI and DiO have fluorescence excitation and emission maxima separated by about 65 nm, facilitating two-color labeling. The emission spectrum of DiA is very broad, allowing it to be detected as green, orange, or even red fluorescence depending on the optical filter used. DiI, DiO, DiD and DiR belong to the dialkylcarbocyanines family of compounds. The spectral properties of the dialkylcarbocyanines are largely independent of the lengths of the alkyl chains, but are instead determined by the heteroatoms in the terminal ring systems and the length of the connecting bridge. They have extremely high extinction coefficients, moderate fluorescence quantum yields, and short excited state lifetimes in lipid environments (~1 ns). The fluorescence spectrum of the dye is shown below.

You can choose the m-Fluoroliposome® based on the type of the fluorescent equipment and filters that you use in your lab. Mannosylated clodronate liposomes cannot be made fluorescent simply due to the potential for inaccurate and/or uninterpretable data being generated by labelled m-Clodrosome®.



Normalized fluorescence emission spectra of DiD, DiI, DiO and DiR (left) and fluorescence excitation and emission spectra of DiA (right).



Macrophage uptake of fluorescent liposome containing DiO.




## FORMULATION INFORMATION

### Five-Color Mannosylated Fluoroliposome® Kit

Lipid Composition	Concentration (mg/ml)	Concentration (mM)	Molar Ratio Percentage
L- $\alpha$ -Phosphatidylcholine	18.8	24.3	70
Cholesterol	4.2	10.9	30
<b>Total</b>	<b>23 mg/ml</b>	<b>35.1 mM</b>	<b>100</b>


Mannosylation	Concentration
4-Aminophenyl- $\alpha$ -D-mannopyranoside	9.53 mol%

Buffer and Liposome Size	Specification
Buffer	Phosphate Buffered Saline
pH	7.4
Liposome Size	1.5-2 $\mu$ m

Fluorescent Dye	Excitation/ Emission (nm)	Concentration (mg/ml)	Concentration (mM)
4-(4-(Dihexadecylamino)styryl)-N-Methylpyridinium Iodide (DiA)	456/590	0.0625	0.0794
			
1,1'-Dioctadecyl-3,3,3',3'-Tetramethylindodicarbocyanine, 4-Chlorobenzenesulfonate Salt (DiD)	644/665	0.0625	0.065
			
1,1'-Dioctadecyl-3,3,3',3'-Tetramethylindodicarbocyanine Perchlorate (DiI)	549/565	0.0625	0.065
			



Fluorescent Dye	Excitation/ Emission (nm)	Concentration (mg/ml)	Concentration (mM)
13,3'-Dilinoleyloxacarboyanine Perchlorate (DiO)	484/501	0.0625	0.071




1,1'-Dioctadecyl-3,3,3',3'-  
Tetramethylindotricarboyanine Iodide (DiR)

750/780

0.00625

0.00651



## TECHNICAL NOTES

- The issue with fluorescent Clodrosome® has to do with the potential for inaccurate and/or uninterpretable data being generated by labelled Clodrosome®. When Clodrosome® induces macrophage apoptosis, the fluorescent lipid incorporated into the Clodrosome® that is disrupted and metabolized in the phagolysosome will be dispersed

among the residual apoptotic bodies which are subsequently phagocytosed by other macrophages. Therefore, fluorescent lipid may be detected in phagocytic cells which never phagocytosed Clodrosome® especially when FACS or fluoroscopy are utilized to detect fluorescent cells (FACS) or fluorescence levels in a tissue homogenate (fluoroscopy). Another potential artifact arises from fluorescent lipid remaining in the extracellular “garbage”, which has not yet been cleared by other phagocytes, generating a high background fluorescence. However, experienced confocal microscopists may be able to differentiate between the punctate fluorescence resulting from fluorescent intact liposomes versus the more diffuse fluorescence characteristic of disrupted liposomes and some have successfully used fluorescent clodronate liposomes to visualize the cellular location of these liposomes by confocal microscopy *in vivo* [2]. A further complicating factor is that published data varies widely as to exactly when clodronate liposomes begin to induce apoptosis in macrophages. Mönkönnen *et al.* show that macrophage death is measurable within the first hour after clodronate liposome treatment on RAW264 cells *in vitro* [3], while many report no signs of macrophage apoptosis until several hours after treatment *in vivo*. The variability in the data is likely due to different liposomal formulations of clodronate as well as the vastly different experimental conditions. Therefore, as with most biological studies, especially those involving liposomes, the amount of time between treating the animal or cells with clodronate liposomes and the onset of apoptosis will need to be established in each experimental model. If the nature of the research demands that Clodrosome® be tracked rather than the control, Encapsula can provide DiI-labelled Clodrosome® upon request, and assuming that the Clodrosome® distribution can definitively be assessed prior to the onset of apoptosis, clear and valid data on the biodistribution of fluorescent Clodrosome® should be obtainable. Still, for most purposes, Fluoroliposome® (fluorescent control liposomes) will provide the required data with far fewer potential artifacts.

- When monitoring monocyte uptake *in vivo* in normal animals, the circulating monocytes may “disappear” or show reduced counts within the first 2 h post-injection due to margination of the monocytes post-liposome phagocytosis. These cells will re-enter the circulation within a few hours. Sunderkötter *et al.* demonstrate this phenomenon and discuss the behavior in detail. Also consider that circulating monocytes have a lifetime of

about 24 h so labeled monocytes will be continually leaving the circulation, even in normal animals, due to aging of the monocytes [4].

- Liposomes may settle when left undisturbed for more than a few hours. Immediately prior to use, in order to ensure a homogeneous liposome suspension, slowly invert the vial several times until the suspension appears homogeneous by visual inspection. Vigorous or erratic shaking will not damage the liposomes, but may induce foaming and bubble formation making it more difficult to accurately measure the desired dosage.
- If the personnel performing intravenous injections are not experienced in or familiar with, precautions for injecting larger volumes (~10% animal weight in ml), viscous liquids or particulate suspensions, consider having extra animals available in case serious injection-related adverse events occur. Dose control animals first to become familiar with large volume injections.
- When dosing intravenously, use standard precautions for dosing larger volumes to animals including the following: a) Warm product to room temperature prior to dosing. b) Ensure that all air bubbles are removed from the syringe prior to dosing; intravenous injection of air bubbles may result in air emboli which can kill or seriously injure animals. c) Inject product at a slow, steady rate of no more than 1 ml/min; decrease infusion rate if animals display any atypical reactions such as unusual agitation.
- Infusion-related adverse reactions usually involve the animal gasping for air or other seizure-like movements. Animals often recover with no apparent permanent injury, but any potential effects on experimental results must be assessed by the researcher.
- Liposomes should be kept at 4 °C and **NEVER** be frozen.

## APPEARANCE

---

m-Fluoroliposome®-DiD, m-Fluoroliposome®-DiI and m-Fluoroliposome®-DiR are blue, pink and dark blue liquid suspensions, respectively. m-Fluoroliposome®-DiA and m-Fluoroliposome®-DiO are yellow liquid suspensions. These Fluoroliposome® products are made of large micron size multilamellar liposomes. Due to their large size, some liposomes



might settle to the bottom of the vial. If left sitting idle in the refrigerator, they will phase separate and form pellets in the bottom of the vial, leaving a clear solution on top. Therefore, the vial should be shaken to form a homogeneous solution prior to use.

## **STORAGE AND SHELF LIFE**

---

### **Storage**

m-Fluoroliposome® products should always be stored at in the dark at 4 °C, except when brought to room temperature for brief periods prior to animal dosing. DO NOT FREEZE. ENS is not responsible for results generated by frozen product.

### **Shelf Life**

m-Fluoroliposome® products are made on daily basis. The batch that is shipped is manufactured on the same day. It is advised to use the products within 60 days of the manufacturing date.

## **REFERENCES AND BACKGROUND READING**

---

1. *Umezawa FA, Eto Y. Liposome targeting to mouse brain: mannose as a recognition marker. Biochemical and biophysical research communications. 1988 Jun 30;153(3):1038-44.*
2. *Polfliet MM, Goede PH, van Kesteren-Hendrikx EM, van Rooijen N, Dijkstra CD, van den Berg TK. A method for the selective depletion of perivascular and meningeal macrophages in the central nervous system. J. Neuroimmunol. 2001 Jun 1;116(2):188–95.*
3. *Mönkkönen J, Liukkonen J, Taskinen M, Heath TD, Urtti A. Studies on liposome formulations for intra-articular delivery of clodronate. Journal of Controlled Release. 1995 Aug;35(2–3):145–54.*

4. [Sunderkötter C, Nikolic T, Dillon MJ, van Rooijen N, Stehling M, Drevets DA, Leenen P. Subpopulations of Mouse Blood Monocytes Differ in Maturation Stage and Inflammatory Response. \*J Immunol.\* 2004 Apr 1;172\(7\):4410–7.](#)
5. [Nagai H, Kuwahira I, Schwenke DO, Tsuchimochi H, Nara A, Ogura S, Sonobe T, Inagaki T, Fujii Y, Yamaguchi R, Wingenfeld L. Pulmonary macrophages attenuate hypoxic pulmonary vasoconstriction via  \$\beta\$ 3AR/iNOS pathway in rats exposed to chronic intermittent hypoxia. \*PLoS One.\* 2015 Jul 1;10\(7\):e0131923.](#)
6. [Zhu Y, Soderblom C, Krishnan V, Ashbaugh J, Bethea JR, Lee JK. Hematogenous macrophage depletion reduces the fibrotic scar and increases axonal growth after spinal cord injury. \*Neurobiology of disease.\* 2015 Feb 28;74:114-25.](#)
7. [Yun MH, Davaapil H, Brockes JP. Recurrent turnover of senescent cells during regeneration of a complex structure. \*Elife.\* 2015;4:e05505.](#)
8. [Arwert EN, Harney AS, Entenberg D, Wang Y, Sahai E, Pollard JW, Condeelis JS. A Unidirectional Transition from Migratory to Perivascular Macrophage Is Required for Tumor Cell Intravasation. \*Cell reports.\* 2018 May 1;23\(5\):1239-48.](#)
9. [Ito T, Ishigami M, Matsushita Y, Hirata M, Matsubara K, Ishikawa T, Hibi H, Ueda M, Hirooka Y, Goto H, Yamamoto A. Secreted Ectodomain of SIGLEC-9 and MCP-1 Synergistically Improve Acute Liver Failure in Rats by Altering Macrophage Polarity. \*Scientific reports.\* 2017 Mar 8;7:44043.](#)
10. [Miron VE, Boyd A, Zhao JW, Yuen TJ, Ruckh JM, Shadrach JL, van Wijngaarden P, Wagers AJ, Williams A, Franklin RJ. M2 microglia and macrophages drive oligodendrocyte differentiation during CNS remyelination. \*Nature neuroscience.\* 2013 Sep;16\(9\):1211.](#)
11. [Andreou K, Sarmiento Soto M, Allen D, Economopoulos V, de Bernardi A, Larkin J, Sibson NR. Anti-Inflammatory Microglia/Macrophages as a Potential Therapeutic Target in Brain Metastasis. \*Frontiers in oncology.\* 2017;7:251.](#)

12. *Alishekevitz D, Gingis-Velitski S, Kaidar-Person O, Gutter-Kapon L, Scherer SD, Raviv Z, Merquiol E, Ben-Nun Y, Miller V, Rachman-Tzemah C, Timaner M. Macrophage-induced lymphangiogenesis and metastasis following paclitaxel chemotherapy is regulated by VEGFR3. Cell reports. 2016 Oct 25;17(5):1344-56.*
13. *Oh SH, Kim HN, Park HJ, Shin JY, Bae EJ, Sunwoo MK, Lee SJ, Lee PH. Mesenchymal stem cells inhibit transmission of  $\alpha$ -synuclein by modulating clathrin-mediated endocytosis in a Parkinsonian model. Cell reports. 2016 Feb 2;14(4):835-49.*
14. *Kano F, Matsubara K, Ueda M, Hibi H, Yamamoto A. Secreted Ectodomain of Sialic Acid-Binding Ig-Like Lectin-9 and Monocyte Chemoattractant Protein-1 Synergistically Regenerate Transected Rat Peripheral Nerves by Altering Macrophage Polarity. STEM CELLS. 2017 Mar 1;35(3):641-53.*